

Control/Tracking Number: 12-SA-3508-ACSM Activity: Scientific Abstract Current Date/Time: 11/1/2011 10:52:06 AM

Validation Of The Bodymetrix Ultrasound System For Percent Body Fat

Author Block: Jean P. Boucher, FACSM¹, Gabriel Chamberland¹, Mylene Aubertin-Leheudre¹, David H. Jones², Ron Rehel², Alain S. Comtois¹. ¹UQAM, Montreal, QC, Canada. ²Concordia University, Montreal, QC, Canada.

Abstract:

2

The Body Metrix (BM) system (IntelaMetrix, CA) is an affordable computer based ultrasound system that provides rapid body composition assessment. We present preleminary data comparing percent body fat (%BF) as measured using the BM system and a DEXA as a reference measurment.

PURPOSE: To establish the validity of the BM computer based ultrasound system for the quantification of %BF. METHODS: %BF was calculated on 10 university aged participants (5 females, 5 males) using the BM system with the Jackson and Pollock 7 sites equation (1978). A DEXA was also used as a reference to measure %BF on the same subjects. The Bland and Altman (1986) statistical procedure was used to determine de validity of the BM %BF.

RESULTS: The level of agreement between both devices was equal to -0.6% BF. As well, no significant correlation was found between the score difference from both devices and the score average of both devices (r = -0.17, p=0.641). The regression analysis between the BM system and the DEXA showed a significant correlation (r=0.71, p=0.023) with a SEE of \pm 3.93 %BF.

CONCLUSIONS: Both devices can be used interchangeably within the range of 10 to 30% BF. The data obtained on a limited number of subjects (n=10) show good agreement and the BM system underestimates by 0.6% on average when compared to the DEXA.

Author Disclosure Information: J.P. Boucher: None.

Category (Complete): 101 body composition Keyword (Complete): Body composition ; Percent body fat Unlabeled/ Investigational Products (Complete): : No Presentation Preference (Complete): Indifferent Area of Interest (Complete): Area of Interest: Applied Science

Additional Info (Complete):

*Do you authorize ACSM to record your presentation?: Yes *Confirmation: I understand and agree to the above terms regarding AV equipment

Payment (Complete): Your credit card order has been processed on Monday 31 October 2011 at 7:54 PM. Status: Complete

OASIS Helpdesk

American College of Sports Medicine 401 West Michigan Street Indianapolis, IN 46202-3233 (317) 637-9200

Leave OASIS Feedback

Powered by <u>OASIS</u>, The Online Abstract Submission and Invitation System SM © 1996 - 2011 <u>Coe-Truman Technologies</u>, Inc. All rights reserved.